


Circuito odontológico Actynox® Aspiração Ativa

Conexão à
espiga da
garrafa



1A Máscara de
Adultos

1B Máscara
Pediátrica

















Montar no
sentido do gás
conforme
indicado















12
Conexão ao
sistema de
aspiração da
cadeira
odontológica



6006156 - Circuito odontológico Actynox® Aspiração Ativa - Adultos - Componentes - Lote: xxxxxx

Item	Cod.	Lot.	Descrição			Quant. Circuito
1A	SIL-MED-1		 CE 0413 Máscara Odontológica Silhouette		Parker Hannifin Corp. 245 Township Line Road, Hatfield, Pennsylvania 19440, United States	5
2	2612001		 CE 0120 Dispositivo standard (Inclui uma bolsa reservatório de 2l e um tubo de condução de oxigénio de 1,8m)		Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK	1
3	1971000		 CE 0120 Conetor 22M - 30F		Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK	1
4	1568000		 CE 0120 Conetor reto 22F - 6MM		Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK	3
5	4222/701		 CE 2797 Filtro bacteriano / vírico		GVS Filter Technology UK Limited NFC House Vickers Industrial Estate Mellishaw Lane Morecambe Lancashire LA3 3EN United Kingdom	5
6	M143		 CE Conetor de Máscara Silhouette da linha inspiratória		Martín Vecino, S.L. c/León Nº80. P.I. Cobo Calleja 28947-Fuenlabrada Madrid – España	1
7	M144		 CE Conetor de Máscara Silhouette da linha expiratória		Martín Vecino, S.L. c/León Nº80. P.I. Cobo Calleja 28947-Fuenlabrada Madrid – España	1
8	MV019		 CE 0123 Linha de condução de oxigénio 1,8 m		SUMBOW MEDICAL INSTRUMENTS CO.,LTD. Add: 14F-01, No.2 Building Shangdong Business Center, #1926 Canghai Road, Yinzhou, 315041 Ningbo. China	2
9	00100510		 CE 0843 Conexão "Y"		QUIROSA, S.A. Crta. Artés a Manresa, km 1,5 P.I. Sta. Mª Artés, Nau 16-17 08271 ARTÉS (Barcelona) Spain	1
10	MV003		 CE 0123 Linha de condução de oxigénio 10 cm		SUMBOW MEDICAL INSTRUMENTS CO.,LTD. Add: 14F-01, No.2 Building Shangdong Business Center, #1926 Canghai Road, Yinzhou, 315041 Ningbo. China	1
11	1950000		 CE 0120 Válvula unidirecional		Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK	1
12	M145		 CE Casquilho de conexão da linha de extração de gases		Martín Vecino, S.L. c/León Nº80. P.I. Cobo Calleja 28947-Fuenlabrada Madrid – España	1

6014259 - Circuito odontológico Actynox® Aspiração Ativa – Pediátrico - Componentes - Lote: xxxxxx

Item	Cod.	Lot.	Descrição			Quant. Circuito
1B	SIL-SM-1		 <p>CE 0413</p> <p>Máscara Odontológica Silhouette</p>		<p>Parker Hannifin Corp. 245 Township Line Road, Hatfield, Pennsylvania 19440, United States</p>	5
2	2612001		 <p>CE 0120</p> <p>Dispositivo standard (Inclui uma bolsa reservatório de 2l e um tubo de condução de oxigénio de 1,8m)</p>		<p>Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK</p>	1
3	1971000		 <p>CE 0120</p> <p>Conetor 22M - 30F</p>		<p>Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK</p>	1
4	1568000		 <p>CE 0120</p> <p>Conetor reto 22F - 6MM</p>		<p>Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK</p>	3
5	4222/701		 <p>CE 2797</p> <p>Filtro bacteriano / vírico</p>		<p>GVS Filter Technology UK Limited NFC House Vickers Industrial Estate Mellishaw Lane Morecambe Lancashire LA3 3EN United Kingdom</p>	5
6	M143		 <p>CE</p> <p>Conetor de Máscara Silhouette da linha inspiratória</p>		<p>Martín Vecino, S.L. c/León Nº80. P.I. Cobo Calleja 28947-Fuenlabrada Madrid – España</p>	1
7	M144		 <p>CE</p> <p>Conetor de Máscara Silhouette da linha expiratória</p>		<p>Martín Vecino, S.L. c/León Nº80. P.I. Cobo Calleja 28947-Fuenlabrada Madrid – España</p>	1
8	MV019		 <p>CE 0123</p> <p>Linha de condução de oxigénio 1,8 m</p>		<p>SUMBOW MEDICAL INSTRUMENTS CO.,LTD. Add: 14F-01, No.2 Building Shangdong Business Center, #1926 Canghai Road, Yinzhou, 315041 Ningbo. China</p>	2
9	00100510		 <p>CE 0843</p> <p>Conexão "Y"</p>		<p>QUIROSA, S.A. Crta. Artés a Manresa, km 1,5 P.I. Sta. Mª Artés, Nau 16-17 08271 ARTÉS (Barcelona) Spain</p>	1
10	MV003		 <p>CE 0123</p> <p>Linha de condução de oxigénio 10 cm</p>		<p>SUMBOW MEDICAL INSTRUMENTS CO.,LTD. Add: 14F-01, No.2 Building Shangdong Business Center, #1926 Canghai Road, Yinzhou, 315041 Ningbo. China</p>	1
11	1950000		 <p>CE 0120</p> <p>Válvula unidirecional</p>		<p>Intersurgical Ltd. Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK</p>	1
12	M145		 <p>CE</p> <p>Casquilho de conexão de linha da extração de gases</p>		<p>Martín Vecino, S.L. c/León Nº80. P.I. Cobo Calleja 28947-Fuenlabrada Madrid – España</p>	1